

Sip/Puff Product Line



Switch Options

Sip/Puff Switch (AC-0304-V2)

The Sip/Puff Switch translates sips and puffs into independent switch closures. The rugged electro-mechanical design requires no external power. Test buttons are integrated to assist in application testing. Electrical interface is provided through a choice of connectors for two standard mono cables or a single stereo cable. A twist-on pneumatic interface provides for direct connection to a range of available user interface options, including Headset and Gooseneck.



Breeze™ (AC-0314)

Breeze is an advanced USB sip/puff switch in an extremely compact package. It connects directly to standard USB ports on computers, laptops, and AAC devices and provides hardware emulation of mouse buttons or joystick functions. For special applications, sensitivity is selectable and dynamic pressure readings can be directly accessed through the USB interface. The integrated pneumatic port connects directly to a range of available user interface options including Headset and Gooseneck.



User Interfaces

Sip/Puff Headset (AC-0305)

Sip/Puff Headset is a comfortable, adjustable head frame. The padded headset slips over the ears and behind the head. This head-referenced design maximizes compatibility with other assistive technologies - like HeadMouse® - that employ head movement to enable a rich interface. A replaceable mouth tube is held in position with an integral stainless steel forming wire. The tube and wire can be bent and cut to adapt to individual needs. A thin, pliable connecting tube interfaces to Sip/Puff Switch or Breeze with a simple twist connector.



Sip/Puff Gooseneck (AC-0315)

Strong, lightweight and versatile, the Sip/Puff Gooseneck features a universal clamp that securely attaches to table tops, beds and wheelchairs. The flexible steel gooseneck and the ball-and-socket mount enable easy and secure mouthpiece positioning. The filtered mouthpiece is replaceable and a thin, pliable connecting tube interfaces to Sip/Puff Switch or Breeze with a simple twist connector.



Sip/Puff Switch Single User Packages

Sip/Puff Switch with Headset (AC-0300-V2)

- Sip/Puff Headset with 6-foot tube
- Sip/Puff Switch
- Cables (1 mono and 1 stereo)
- 4 mouth tubes

Sip/Puff Breeze with Headset (AC-0313-H)

- Sip/Puff Headset with 6-foot tube
- Breeze
- 4 mouth tubes

Sip/Puff Switch with Gooseneck (AC-0300-G)

- Sip/Puff Gooseneck with clamp
- Sip/Puff Switch
- 6-foot tube
- Cables (1 mono and 1 stereo)
- 4 mouth pieces

Sip/Puff Breeze with Gooseneck (AC-0313-G)

- Sip/Puff Gooseneck with clamp
- Breeze
- 6-foot tube
- 4 mouth pieces

Sip/Puff Switch Multi-User Package

Multi-user Sip/Puff Switch (AC-0300-MUV2)

The multi-user version of the Sip/Puff Switch provides healthcare institutions, assessment centers, hospitals, etc. with a means to support multiple users. The comfortable, adjustable headset affords comfort and convenience. The disposable mouthpiece includes an in-line aseptic filter to enable hygienic, worry-free use by a succession of users. The Multi-user Sip/Puff Switch is packaged in a compact case for storage and transport.



Origin Instruments Corporation 854 Greenview Drive GrandPrairie, TX 75050
972.606.8740 (Phone) 972.606.8741(Fax)

www.orin.com Shop online at: <http://shop.orin.com>