

## Notebook Bracket Installation Instructions

The Notebook Bracket is for use with notebook computers, some desktop computers with flat-panel displays (especially where the display frame is too narrow for fastening Dual Lock), and some augmentative communication devices (AAC).

The bracket comes with low-density Dual Lock and rubber pads preinstalled for securely attaching the HeadMouse Extreme to the bracket. The HME is positioned with its high-density Dual Lock fastener aligned with the Dual Lock fastener on the bracket. The HME should be oriented such that the connectors point down. Firmly press the HME against the rubber pads, then press back until the Dual Lock mates. There may or may not be a clicking sound as the Dual Lock fasteners connect.



Notebook Bracket



Bracket Mounted on Notebook Computer

Plug the small end of the straight USB cable supplied with your bracket into the corresponding connector on the HeadMouse Extreme.

The bracket comes with two sets of wire clips. The shorter clips are used with narrow displays, while the longer clips are appropriate for a thicker display monitor. Loosen the thumbscrews on either side of the bracket, and slide out the wire clips far enough to fit over the top edge of your computer display or AAC device. If the clips are not long enough, replace with the longer clips. Once loosely on the display, push gently on the clip to firmly grip the display (though not so tight as to warp your display), while you tighten the thumbscrews.



Front View on Notebook Computer